


PHYSIOTHERAPEUTIC COLD THERAPY AND HOT THERAPY



Analgesic and anti-inflammatory –  
relaxing and circulation-promoting





**Ice Compresses for cold therapy**

- Safety ensured by double layers of plastic
- Special gel facilitates accurately dosed
- Cold emission
- Flexible up to -20° C
- Easy to use: store in freezer
- Hygienic, can be disinfected
- Durable and re-usable (cost-efficient)

Thermoform® Ice Compresses were developed for treatment in cases in which cold therapy is indicated. The durable, re-usable compress is filled with a **special gel** that stores cold for an **extremely long period** and ensures a precisely dosed emission of cold. In addition, the Thermoform® Ice Compresses remain **flexible up to -20° C**, which makes handling easier and facilitates optimal adaptation to the body parts to be treated. Simply place it in the freezer to cool it down.

**Indications include:**

- Muscle and joint inflammations
- Rheumatic disorders
- Bruises, Sprains, Pulled muscles
- Fresh injuries
- Nosebleeds
- Vascular inflammations
- Thromboses
- Headaches/migraines

**Contraindications:**

- Impaired arterial blood flow
- Raynaud's syndrome (blue fingertips)



**Thermoform® # 1110**  
Ice Compress  
Size: 22 x 22 cm

**Thermoform® # 1120**  
Ice Compress  
Size: 24 x 32 cm

**Thermoform® # 1130**  
Ice Compress  
Size: 24 x 42 cm

**Thermoform®**  
Ice Collar with 2 tapes for the neck and joints,  
· for adults # 1140  
Size: 10 x 30 cm  
· for children # 1150  
Size: 7 x 25 cm

**Thermoform® # 1190**  
Shoulder Ice Compress with tapes

**Thermoform® # 1120 and 1130**  
application example:  
for cooling elbows at Epicondylitiy, ankle etc. (Fig. with fixation bandage)



# COLD THERAPY PRODUCTS FOR CLINICS AND PRIVATE PRACTICE



## IceCap to prevent hair loss during chemotherapy

- Unisize
- Isolation hood for constant cold
- Cold compress is interchangeable

In many cases chemotherapies used in the treatment of malignant tumors and metastases have an adverse side effect: patients lose all their hair. For women in particular, this causes psychological stress above and beyond that of the illness itself.

Targeted hypothermia of the scalp and the resulting ischemia, caused by the con-


traction of the capillaries supplying the scalp with blood, reduces the effect of the medications applied in this area to a minimum. In most cases, hypothermic treatment with the Thermoform® IceCap is successful.

The Thermoform® IceCaps are ideally adapted to the patient's **anatomical features** and come in all head sizes.

**Thermoform® # 1180**  
IceCap Set with Velcro tape and insulation hood

**Thermoform® # 1181**  
IceCap with Velcro tape

**Thermoform® # 1182**  
IceCap Insulation hood (separate)



**Thermoform® # 1190 and 1140**  
application example:

for cooling of the shoulder joint  
Capsulitis, arthritis, bursitis etc.  
and for cooling in throat inflammation  
or after tonsillectomy



## Paraffin Compresses for heat therapy

- High wear resistance
- Hygienic (easy to clean)
- High resistance to oils and solvents
- Easy to use

Thermoform® Paraffin Compresses facilitate practical, hygienic and affordable heat therapy in professional practice or in the patient's home.

A high-quality, **extremely flexible plastic material** is used in all Thermoform® Paraffin Compresses. This allows a safe, convenient and comfortable use of the products. Paraffin has a unique, outstanding heating effect. The specific melting point of paraffin ensures **particularly**

**consistent heat dissemination** throughout the entire treatment period.

**Indications include:**

- Circulatory disorders
- Bruises, pulled muscles
- Sprains
- Spinal and intervertebral disk damage
- Tense muscles
- Sciatica
- Rheumatic illnesses



**Thermoform® # 1020**  
Paraffin Compress  
Size: 24 x 32 cm


**Thermoform® # 1030**  
Paraffin Compress  
Size: 24 x 42 cm

**Thermoform® # 1040**  
Paraffin Compress  
with 2 tapes for the  
neck and joints  
Size: 10 x 30 cm

**Thermoform® # 1050**  
Paraffin Compress  
with 4 tapes for the  
shoulder area



## Moor- and Fango Compress for traditional pain relief



Moor mud is a traditional ingredient in therapeutic heat treatments. Thanks to its high waterbinding and heat storage capacities, the heated moor mud can be tolerated on the skin at higher temperatures than water, and emits its warmth over a longer period of time. We recommend article # 1030 as a heat transfer medium.

**Thermoform® # 1220**  
**Moor Compress (disposable)**  
Rich in humic acids, estrogens, tannic acid and mineral salts. Apply as partial moor pack directly to the skin.  
Size: approx. 30 x 40 cm

**Thermoform® # 1221**  
**Fango Compress (disposable)**  
Volcanic mineral mud for partial fango packs for the treatment of rheumatic diseases or disorders of the inner organs.  
Size: approx. 40 x 60 cm



## Warm/Cold Compresses for pain relief

- Durable nylon fabric
- Simple to handle
- Double funktion

The cold / warm compress made of high-quality, tough nylon fabric. Durable, re-usable, easy to care for. To use as a cold compress, cool in freezer (up to -20 °C). To use as a warm compress, heat in hot water or microwave.

### Thermoform® Softy-Pack

- Size: 12 x 15 cm # 1000
- Size: 15 x 24 cm # 1001
- Size: 24 x 30 cm # 1002



## Hot water bottle with special gel for the microwave

- Comfortable warmth in just a few minutes
- Relaxes the muscles
- Pleasant film surface

Wherever warmth contributes to well-being, these cushions are indispensable aids. In just a few minutes, cozy, pleasant warmth is created quickly and easily. Foil colour: red

### Thermo-Star 1

#### Microwaveable hot water bottles

- round, Ø = 25 cm # 1502
- Matching cover, blue # 1505
- Hot water bottle + cover # 1500



## Thermoform® accessories



### Thermoform® Covers (white)

- for Ice Compress 1110 # 1401
- for Ice Compress 1120, 1160 # 1402
- for Ice Compress 1130 # 1403
- for Ice Compress 1140 # 1404

### Thermoform®

#### Retention bandages for ice compresses

- 90 x 7,5 cm # 1303
- 120 x 15 cm # 1304
- 100 x 11 cm # 1305 (with pocket)




# 1305

## HEADQUARTER

---

You will find us in Garbsen, close to Hannover (Northern Germany) with very good access to the freeway and the airport Hannover!



**FACEBOOK**  
B.u.W. Schmidt GmbH




**YOUTUBE**  
B.u.W. Schmidt



**DOWNLOAD**  
[www.buwschmidt.com](http://www.buwschmidt.com)

---

Visit our website for further information and product brochures for downloading.



Since 1961, B. u. W. Schmidt GmbH in Garbsen near Hanover has produced precision products made of high-frequency welded flat sheets for the health and emergency services sector. Faultless quality and workmanship, easy handling and rugged materials are the basis for the uncompromising safety of our products.

Our service ensures short delivery times and fast repairs in our factory. The quality management system of B.u.W. Schmidt GmbH is EN ISO certified. All equipment complies with the German Medical Device Act (MPG) and bears the CE mark.

If you cannot find the desired items or dimensions in this brochure, please contact us regarding custom production. Our computer-aided production system enables us to make equipment to meet your individual requirements in all sizes up to 150 x 250 cm.



**Schmidt**

PLASTIC PRODUCTS  
for health and emergency services

---

B. u. W. Schmidt GmbH ■ Porschestraße 29 ■ D-30827 Garbsen  
Phone: +49 (0)5131/6008 ■ Fax: +49 (0)5131/96393 ■ E-Mail: [mail@buwschmidt.com](mailto:mail@buwschmidt.com) ■ [www.buwschmidt.com](http://www.buwschmidt.com)